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REF

560015

10 Röhrrchen

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<http://www.medco.eu/pdf/gbs-ablesehilfe.pdf>

Literatur

- Schrag S, Gorwitz R, Fultz-Butts K, Schuchat A (2002): Prevention of perinatal group B streptococcal disease. Revised guidelines from CDC. MWR Recomm Rep 51: 1-22.
- Schuchat A: Group B streptococcus. Lancet 1999, 353:51.
- Kolben M et al.: Peripartales Management bei mütterlicher Streptokokken B-Kolonisation. Pädiat Prax 1998,54:21-26.
- Boyer KM et al.: Prevention of early-onset neonatal group B streptococcal disease with selective intrapartum chemoprophylaxis. NEJM 1986,314:1665-1669.
- Milatovic D et al.: B-Streptokokken-Screening mittels GBS-Medium in der Geburtshilfe. Immun Infekt 1995,23:134-136.
- Hafner E et al.: Group B streptococci during pregnancy: a comparison of two screening and treatment protocols. Am J Obstet Gynecol 1998, 178:667-681.
- Islam AKM: Rapid recognition of group B streptococci. Lancet 1977, 256-7.
- Schuchat A: Epidemiology of group B streptococcal disease in the United States: shifting paradigms. Clin Microbiol Rev 1998, 11:497-513.
- Votava et al.: Use of GBS media for rapid detection of Group B streptococci in vaginal and rectal swabs from women in labor. Eur J Clin Microbiol Infect Dis 2001; 20:120-12.

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Bestell-Nr.: 560015

GBS – Röhrcchen 10 Stück