COMPACTDRY BC

Instructions

Sample Preparation

Viable count in water or liquid foodstuff

1 mL of sample (diluted if necessary) is deposited on the BC CompactDry plate.

Viable count in solid food

Add buffer solution to the sample and homogenize using a Stomacher® (as per manufacturer instructions). 1 mL of sample (diluted if necessary) is dispensed in the centre on the BC CompactDry plate.

Viable count in swab test sample

Use the swab to wipe the surface, put the swab back into the tube and screw tightly. Shake the tube to mix. 1 mL of sample (diluted if necessary) is dispensed in the centre on the BC CompactDry plate.

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
CompactDry Protocol

1 Open the cap **2** Dispense 1 mL of sample on the middle of the CompactDry plate **3** Sample diffuses automatically and evenly across the sheet transforming the dry medium into a gel within seconds. **4** Replace the cap, label and turn the plate over and Incubate for the appropriate time and **5** required temperature. **6** Following incubation, count the number of coloured colonies.

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Incubation temperature 30 +/- 1 °C

Incubation time:

First inspection: 24 +/- 2 hours plus a further 24 +/- 2 hours

Interpretation

Bacillus cereus forms light blue to blue colonies.

Other potential contaminated bacteria may also grow and form white colonies in this plate. However, only light blue/blue colonies should be counted. Related species, Bacillus thuringiensis may also grow and form light blue/blue colonies.

General information

The growth area is 20 cm². The back of the plate has a 1 cm² grid to ease colony counting. If there are a large number of colonies, the total viable count can be obtained by multiplying the average number of colonies in several grids by 20.

High concentrations (> 300 CFU) will cause the entire growth area to become blue. In this case dilute the sample.

After use please follow local disposal regulations.

Storage and shelf life

Store at room temperature (+1 to +30 °C). Shelf life 18 months after manufacturing.

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