Storage up  
to 30°C

18 months

## Sample Preparation

### Viable count in water or liquid foodstuff

1 mL of sample (diluted if necessary) is deposited on the PA CompactDry plate.

### Viable count using filter technology

See filter method instructions ([www.nissui-ps.com](http://www.nissui-ps.com)).

### Viable count in solid foodstuff


Add appropriate volume of buffer solution (e.g. PBS, 1 part solid foodstuff + 9 parts buffer) to the sample and homogenize using a stomacher®. 1 mL of sample (diluted if necessary) will be dispensed in the middle on the PA CompactDry plate.

### Viable count in swab test specimen

Use the swab to wipe the surface, put the swab back into the tube and screw tightly. Shake the tube to mix. 1 mL of sample (diluted if necessary) will be dispensed in the middle on the PA CompactDry plate.

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


## CompactDry Protocol

**1** Open the cap **2** Dispense 1 mL of sample on the middle of the CompactDry plate or add filter as instructed **3** Sample diffuses automatically and evenly across the sheet transforming the dry medium into a gel within seconds. **4** Replace the cap, label and turn the plate over. **5** Incubate for the required time and temperature. **6** Following incubation, count the number of coloured colonies.

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**Incubation temperature 34-36 °C**  
**Incubation time 45 - 51 hours**

#### **Interpretation**

P. aeruginosa forms red colonies which are generally surrounded by a greenish-yellow halo. Tests have shown that around 20% of P. aeruginosa strains form no greenish-yellow halo.

#### **General information**

Some Pseudomonas spp. and Alcaligenes spp. may also create red colonies but are never surrounded by a greenish-yellow halo.

Antibiotics in the nutrient medium of CompactDry PA will inhibit the growth of bacteria others than the species named above.

The growth area is 20 cm<sup>2</sup>. The back of the plate has a 1 cm<sup>2</sup> grid to ease colony counting. If there are a large number of colonies, the total viable count can be obtained by multiplying the average number of colonies in several grids by 20.

High concentrations (>200 CFU) will cause the entire growth area to become red or greenish-yellow. In this case dilute the sample.

After use please follow local disposal regulations.

#### **Storage and shelf life**

Keep at room temperature (+ 1 - + 30°C). Shelf life 19 months after manufacturing.

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