# COMPACTDRY YM

## Instructions

### Sample Preparation

#### Viable count in water or liquid foodstuff

1 mL of sample (diluted if necessary) is deposited on the YM CompactDry plate.

#### Viable count in solid food

Add buffer solution to the sample and homogenize by Stomacher® (as per manufacturer instructions). 1 mL of sample (diluted if necessary) is dispensed in the centre on the YM CompactDry plate.

#### Viable count in swab test specimen

Use the swab to wipe the surface, put the swab back into the tube and screw tightly. Shake the tube to mix. 1 mL of sample (diluted if necessary) is dispensed in the centre on the YM CompactDry plate.

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
## CompactDry Protocol

**1** Open the cap **2** Dispense 1 mL of sample on the centre of the CompactDry plate. **3** Sample diffuses automatically and evenly across the sheet transforming the dry sheet into a gel within seconds. **4** Replace the cap, label and turn the plate over and **5** Incubate for the appropriate time and at required temperature. **6** Following incubation, count the number of coloured colonies.

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Incubation temperature 25 - 30 °C  
Incubation time 3 - 7 days

**Interpretation**

Blue coloured colonies indicate yeasts. Mould form cottony colonies with characteristic colours.

**General information**

Certain yeasts do not form blue colonies.

Antibiotics in the medium inhibit bacterial growth.

The growth area is 20 cm<sup>2</sup>. The back of the plate has a 1 cm<sup>2</sup> grid to ease colony counting. If there are a large number of colonies, the total viable count can be obtained by multiplying the average number of colonies in several grids by 20.

High concentrations (> 300 CFU) will cause the entire growth area to become blue/green. In this case dilute the sample.

After use please follow local disposal regulations.

**Storage and shelf life**

Store at room temperature (+1 - +30 °C). Shelf life 18 months after manufacturing.

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