
COMPACTDRY YMR

Instructions

Sample Preparation

Viable count in water or liquid foodstuff

1 mL of sample (diluted if necessary) is deposited on the YMR CompactDry plate.

Viable count in solid food

Add 9 part of buffer solution (0.1 % peptone water or MRD) to 1 part of food sample and homogenize by stomacher® (as per manufacturer instructions). 1 mL of sample (diluted if necessary) is dispensed in the centre on the YMR CompactDry plate.

Viable count in swab test specimen

Use the swab to wipe the surface, put the swab back into the tube and screw tightly. Shake the tube to mix. 1 mL of sample (diluted if necessary) is dispensed in the centre on the YMR CompactDry plate.

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
CompactDry Protocol

1 Open the cap **2** Dispense 1 mL of sample on the centre of the CompactDry plate. **3** Sample diffuses automatically and evenly across the sheet transforming the dry sheet into a gel within seconds. **4** Replace the cap, label and turn the plate over and **5** Incubate for the appropriate time and at required temperature. **6** Following incubation, count the number of coloured colonies.

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Incubation temperature 25 - 30 °C
Incubation time 48 - 72 hours

Interpretation

Blue coloured colonies indicate yeasts. Mould form cottony colonies with characteristic colours.

General information

Certain yeasts do not form blue colonies.

Antibiotics in the medium inhibit bacterial growth.

The growth area is 20 cm². The back of the plate has a 1 cm² grid to ease colony counting. If there are a large number of colonies, the total viable count can be obtained by multiplying the average number of colonies in several grids by 20.

High concentrations (> 300 CFU) will cause the entire growth area to become green/blue. In this case dilute the sample.

After use please follow local disposal regulations.

Storage and shelf life

Store at room temperature (+1 - +30 °C).

Shelf life 18 months after manufacturing.

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NISSUI PHARMACEUTICAL CO.,LTD. (Japan)